


• Indication:	Titanbasis 2. Generation (local-milling)
• DME Version:	MEDENTIKA_Ti-Base_2.Generation_all_series_local-milling_V20-10_18-20.dme
• 3S-DentalSystem:	2018-2020

- **Possibility of rotating of the Ti-Base**

You can rotate the Ti-Base (all series) around the own implant axis. With this option you don't need to rescan a case, when the Scanbody was inserted in a wrong position. The direction of the Scanbody (flatside) is no longer important like in the versions before.



- **New Scanbody components**

The implementation of the Elos Accurate® Scanbodies and 3Shape Scanbodies gives you more options to chose within the portfolio of MEDENTIKA®.

