

Kompatibilitätsname	Serie	Interface	Medentika Scanbody	ELOS-Scanbody IO	ELOS-Scanbody IO SA	ELOS-Scanbody
Bredent Medical / SKY**	B-Serie	RP/NP	B 1400			
BEGO Implant Systems/ Semados** S/SC/SCX/RS/RSX/RI	BS-Serie	D 3,25-5,5	BS 1400			
Altatec / Camlog**	C-Serie	D 3,3 D 3,8 D 4,3 D 5,0 D 6,0	C 1400 C 1410 C 1420 C 1430 C 1430	IO 9A-A IO 9A-B IO 9A-C IO 9A-D IO 9A-D		9A-A 9A-B 9A-C 9A-D 9A-D
Medentis Medical / ICX	CX-Serie	D 3,75-4,8	CX 1400	IO_20A-A		20A-A
Altatec / Conelog**	D-Serie	D 3,3 D 3,8-4,3	D 1400 D 1410			
Dentium / SuperLine, Implantium, Implantium II	DT-Serie	D 3,6-7,0	DT 1400-L			
Nobel Biocare / NobelReplace** Tapered	E-Serie	NP 3,5 RP 4,3 WP 5,0 D 6,0	E 1400 E 1410 E 1420 E 1430	IO 2A-A IO 2A-B IO 2A-C IO 2A-D		2A-A 2A-B 2A-C 2A-D
DENTSPLY Implants / ASTRA TECH OsseoSpeed** EV	EV-Serie	D 3,0 D 3,6 D 4,2 D 4,8 D 5,4	EV 1400 EV 1410 EV 1420 EV 1430 EV 1440	IO 3B-A IO 3B-B IO 3B-C IO 3B-D IO 3B-E		
Nobel Biocare/NobelActive** NobelReplace**Conical	F-Serie	D 3,0 NP 3,5 RP 4,3-5,0	F 1420 F 1400 F 1410	IO 2B-D IO 2B-A IO 2B-B	IO 2B-D SA IO 2B-A SA IO 2B-B SA	2B-D 2B-A 2B-B
BIOMET 3i / Certain**	H-Serie	D 3,4 D 4,1-5,0	H 1400 H 1410	IO 7B-A IO 7B-B		7B-A 7B-B
BIOMET 3i / External Hex	I-Serie	D 3,4 D 4,1-5,0	I 1400 I 1410	IO 7A-A IO 7A-B		7A-A 7A-B
Nobel Biocare / Brånemark System**	K-Serie	NP 3,5 RP 4,1 WP 5,1	K 1400 K 1410 K 1420	IO 6A-A IO 6A-B IO 6A-C		6A-A 6A-B 6A-C
Straumann / Bone Level	L-Serie	NC 3,3 RC 4,1-4,8	L 1400 L 1410	IO 4A-A IO 4A-B		4A-A 4A-B
Straumann / Tissue Level	N-Serie	NNC 3,5 RN 4,8 WN 6,5	N 1430 N 1410 N 1420	IO 4C-A IO 4B-A IO 4B-B		4C-A 4B-A 4B-B
NEOSS	NE-SERIE	SP D3,5-6,5	NE 1400-L	IO 1A-A		1A-A
OSSTEM Implants TS System	OT-Serie	Mini Regular	OT 1400-L OT 1410-L			
Zimmer Dental Tapered Screw-Vent** MIS SEVEN Internal Hex	R-Serie	D 3,5 D 4,5 D 5,7	R 1400 R 1410 R 1420	IO 11A-B IO 11A-C IO 11A-D		11A-B 11A-C 11A-D
DENTSPLY SIRONA ASTRA TECH OsseoSpeed** TX	S-Serie	D 3,0 D 3,5-4,0 D 4,5-5,0	S 1430 S 1400 S 1420	IO 3A-A IO 3A-B IO 3A-C		3A-A 3A-B 3A-C
DENTSPLY SIRONA XIVE** S	T-Serie	D 3,4 D 3,8 D 4,5-5,5	T 1400 T 1405 T 1410	IO 8A-B IO 8A-C IO 8A-D		8A-B 8A-C 8A-D
DENTSPLY Sirona ANKYLOS** C/X	Y-Serie	C/X 3,5-7,0	Y 1400X-L	IO 8B-A		8B-A
Multi-unit	all series	NP/RP MUA45	8500/8513	IO 2C-A		2C-A