

Indication:	Multi-unit (local-milling)
Version:	MEDENTiKA_Multi-unit_all-series_local-milling_V20-09.zip
exo-DentalSystem:	all systems

- **New Titanium Base Flex: TiBase 0-31-20/8519**

Flexible chimmney heights 4,5mm – 14,5mm (cut able in 2mm-steps).



- **Titanium Base ASC: 6-31-02 (ANGULATED SCREWCHANNEL)**

The screw channel of the TiBase 6-31-02 in combination with the Ball Torx screw 6-31-01, is circular angulated up to 25°.

